EluNIR-PERL™

Ridaforolimus Eluting Coronary Stent System

[EN] INSTRUCTIONS FOR USE

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[ES] INSTRUCCIONES DE USO

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TABLE OF CONTENTS

3	SYMBOLS USED IN LABELING		
	1 PRODUCT DESCRIPTION		
	1.1 DEVICE COMPONENTS DESCRIPTION		
	1.2 DRUG COMPONENT DESCRIPTION		
	1.3 EXPECTED LIFETIME OF THE DEVICE		
	2 INTENDED PURPOSE AND INDICATIONS FOR USE		
	2.1 INTENDED PURPOSE		
	2.2 INDICATIONS FOR USE		
4	3 CONTRAINDICATIONS		
	4 WARNINGS		
	5 PRECAUTIONS		
	5.1 GENERAL PRECAUTIONS		
	5.2 ANTIPLATELET REGIMEN		
	5.3 USE OF MULTIPLE STENTS		
	5.4 USE IN CONJUNCTION WITH OTHER PROCEDURES		
	5.5 BRACHYTHERAPY		
	5.6 USE IN SPECIAL POPULATIONS		
	5.7 LESION/VESSEL CHARACTERISTICS		
	5.8 DRUG INTERACTIONS		
5	5.9 MAGNETIC RESONANCE IMAGING (MRI) SAFETY INFORMATION		
	5.10 STENT HANDLING		
	5.11 STENT PLACEMENT		
	5.12 STENT SYSTEM REMOVAL		
	5.13 POST PROCEDURE		
	6 DRUG INFORMATION	6	
	6.1 MECHANISM OF ACTION		
	6.2 DRUG INTERACTIONS		
	6.3 MUTAGENESIS, CARCINOGENICITY AND REPRODUCTIVE TOXICITY		
	7 POTENTIAL ADVERSE EVENTS/COMPLICATIONS		
	8 INDIVIDUALIZATION OF TREATMENT		
	8.1 PATIENT COUNSELING AND PATIENT INFORMATION		
	8.2 CLINICAL BENEFITS		
	9 HOW SUPPLIED		
	10 OPERATOR'S MANUAL	7	
	10.1 INSPECTION PRIOR TO USE		
	10.2 MATERIALS REQUIRED		
	10.3 INCIDENTS REPORT		
	10.4 PREPARATION		
	10.5 DELIVERY PROCEDURE		
	10.6 DEPLOYMENT PROCEDURE		
	10.7 REMOVAL PROCEDURE		
	10.8 POST-DEPLOYMENT DILATATION OF STENT SEGMENTS		
	10.9 PATIENT IMPLANT CARD & PATIENT INFORMATION GUIDE		
	11 REUSE PRECAUTION STATEMENT	8	
	12 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE		
	13 DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY	9	

EXPLANATION OF SYMBOLS ON LABELS AND PACKAGING

	Manufacturer		Min. Guiding Catheter ID
	Authorized Representative in the European Community		Diameter
	Consult Instructions for Use		Balloon Pressure
	Do Not Use if Package is Damaged		Nominal Pressure
	Keep Dry		Rated Burst Pressure
	Keep Away from Sunlight		Do Not Exceed Rated Burst Pressure
	Sterilized Using Ethylene Oxide		Catalogue Number
	Do Not Reuse		Lot Number
	Use By		Medical Device
	MR Conditional		Unique Device Identifier
	Do Not Resterilize		25°C / 77°F with transient excursions allowed up to 40°C / 104°F
	Medicinal Substance		Single sterile barrier system with two protective packaging outside
	Contains Hazardous Substances		Single sterile barrier system with protective packaging outside
	Nonpyrogenic		Single sterile barrier system

INSTRUCTIONS FOR USE

EluNIR-PERL™ Ridaforolimus Eluting Coronary Stent System

A paper copy of the IFU can be requested from contact@medinol.com. The paper copy will be delivered, free of charge, within seven days.

EluNIR-PERL™ Ridaforolimus Eluting Coronary Stent System (EluNIR-PERL™) is provided sterile for single use only.

1 PRODUCT DESCRIPTION

The Medinol EluNIR-PERL™ Ridaforolimus Eluting Coronary Stent System (EluNIR-PERL™) is a single use device/ drug combination product composed of the following device components: the coronary stent and its delivery system, and a drug component (a formulation of ridaforolimus in a polymer coating blend). The EluNIR-PERL™ characteristics are described in **Table 1**.

TABLE 1: EluNIR-PERL™ RIDAFOROLIMUS ELUTING CORONARY STENT SYSTEM PRODUCT DESCRIPTION

AVAILABLE STENT LENGTHS (mm)	8, 12, 15, 17, 20, 24, 28, 33, 38, 44	
AVAILABLE STENT DIAMETERS (mm)	2.25*, 2.5*, 2.75, 3.0, 3.5, 4.0	
STENT MATERIAL	A medical grade L-605 Cobalt Chromium [CoCr], annealed, with platinum markers.	
DRUG COMPONENT	A coating of polymers loaded with ridaforolimus in a formulation applied to the entire surface of the stent at a dose of approximately 1.1 µg/mm²	
DELIVERY SYSTEM WORKING LENGTH (cm)	140 cm	
DELIVERY SYSTEM DESIGN	Single access port to inflation lumen; guidewire exit notch (RX-Port) is located 30 cm from distal tip; designed for guidewires ≤ 0.014" [0.36 mm]	
STENT DELIVERY SYSTEM	Expandable balloon with two radiopaque markers to indicate balloon positioning and expanded stent length	
BALLOON INFLATION PRESSURE	Nominal pressure: For diameter 2.25 mm: 8 atm [811 kPa] For diameters 2.5-4.0 mm: 10 atm [1013 kPa] Rated Burst Pressure (RBP), for all diameters: 18 atm [1824 kPa]	
MINIMUM GUIDING CATHETER INNER DIAMETER	≥5F [0.056"/ 1.42 mm]	
CATHETER SHAFT OUTER DIAMETER	Proximal	2.1F [0.69mm]
	Distal	2.7F [0.90 mm] for products of 8 mm - 28 mm length
		2.9F [0.97 mm] for products of 33 mm - 44 mm length

* The 2.25 & 2.5 mm diameter stents for the EluNIR-PERL™ stent system are available in lengths up to 33 mm long.

1.1 DEVICE COMPONENTS DESCRIPTION

The EluNIR-PERL™ stent system consists of a ridaforolimus eluting coronary stent component premounted onto an RX delivery system. The stents are made from a cobalt-based alloy and are coated with a drug/ polymer coating, which consists of a Poly n-Butyl Methacrylate (PBMA) polymer, a CarboSil®20 55D polymer and the active pharmaceutical ingredient (API) ridaforolimus. The EluNIR-PERL™ stents have a drug density of 1.1 µg/mm². The stent has two radio-opaque platinum markers at either end to enhance radiopacity. The EluNIR-PERL™ delivery system provides a means of delivering the stent through the coronary vasculature and, once in the desired location, expands the stent through balloon inflation. The catheter has a hydrophilic coating on the outer surface of the distal shaft.

1.2 DRUG COMPONENT DESCRIPTION

The coating on the EluNIR-PERL™ stent consists of a polymer coating blend [Poly n-Butyl Methacrylate (PBMA) and CarboSil® 20 55D, (inactive ingredients)], and the active pharmaceutical ingredient (API) ridaforolimus. These are the same active pharmaceutical ingredient and polymers as are used in the EluNIR™ Ridaforolimus Eluting Coronary Stent System [EluNIR™].

1.3 EXPECTED LIFETIME OF THE DEVICE

The expected lifetime of the EluNIR-PERL Ridaforolimus Eluting Coronary Stent System is a minimum of 10 years. EluNIR-PERL Ridaforolimus Eluting Coronary Stents are not intended to be removed after implantation. The expected lifetime was established based on extensive DV testing performed according to state-of the art standards including durability and fatigue testing, long-term animal studies and clinical studies with long-term follow-up of 5 years.

2 INTENDED PURPOSE AND INDICATIONS FOR USE

2.1 INTENDED PURPOSE

The EluNIR-PERL™ Ridaforolimus Eluting Coronary Stent System is a single use device/ drug combination product, intended to improve coronary artery luminal diameter in patients with symptomatic heart disease due to de novo lesions. The EluNIR-PERL™ delivery system provides a means of delivering the stent through the coronary vasculature and, once in the desired location, expands the stent through balloon inflation.

2.2 INDICATIONS FOR USE

The EluNIR-PERL™ Ridaforolimus Eluting Coronary Stent System is indicated for improving coronary luminal diameter in patients with symptomatic heart disease due to de novo lesions of length ≤42 mm in native coronary arteries with reference diameters of 2.25mm to 4.25mm and according to ACC/AHA therapeutic guidelines.

3 CONTRAINDICATIONS

Coronary artery stenting is generally contraindicated in the following patient types:

- Patients who cannot receive recommended antiplatelet and/or anticoagulation therapy
- Patients judged to have a lesion which prevents complete inflation of an angioplasty balloon or proper placement of the stent or delivery system.
- Patients with hypersensitivity or allergies to aspirin, heparin, clopidogrel, ticlopidine, drugs such as ridaforolimus or similar drugs, or any other analogue or derivative, polymers, cobalt, chromium, nickel, molybdenum, or contrast media.

4 WARNINGS

- Please ensure that the inner package has not been opened or damaged as this would indicate that the sterile barrier has been breached.
- The use of this device carries the associated risks of thrombosis, vascular complications and/or bleeding events.
- This product should not be used in patients who are not likely to comply with the recommended antiplatelet therapy.

5 PRECAUTIONS

5.1 GENERAL PRECAUTIONS

- Stent implantation should only be performed by interventional cardiologists who have been trained in cardiovascular interventional techniques, in a cardiac catheterization laboratory or a surgical suite.
- Subsequent restenosis may require repeat dilatation of the arterial segment containing the stent. Long-term outcomes following repeat dilatation of the stent is presently not well characterized.
- Risks and benefits should be considered in patients with severe reaction to contrast agent.
- Patients with known hypersensitivity to the product components (stent substrate, polymer[s], drug substance) may suffer an allergic reaction to this implant.
- Do not expose or wipe the product with organic solvents such as alcohol.
- Care should be taken to control the guiding catheter tip during stent delivery, deployment, and balloon withdrawal. Before withdrawing the stent delivery system, visually confirm complete balloon deflation by fluoroscopy to avoid guiding catheter movement into the vessel and subsequent arterial damage.
- Stent thrombosis is a low-frequency event that is frequently associated with myocardial infarction (MI) or death.
- When DES are used outside the specified Indications for Use, patient outcomes may differ from the results observed in the EluNIR™ clinical trials.
- Compared to use within the specified Indications for Use, the use of DES in patients and lesions outside of the labeled indications, including more tortuous anatomy, may have an increased risk of adverse events, including stent thrombosis, stent embolization, MI, or death.
- One or more components of this device contains the following substance defined as CMR 1B in a concentration >0.1% weight/weight.: Cobalt (CAS No. 7440-48-4, EC No. 231-158-0)
Current scientific evidence supports that medical devices manufactured from cobalt alloys and stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

5.2 ANTIPLATELET REGIMEN

Dual antiplatelet therapy (DAPT) using a combination treatment of aspirin with a P2Y12 platelet inhibitor after percutaneous coronary intervention (PCI), has been shown to reduce the risk of stent thrombosis and ischemic cardiac events, but increases the risk of bleeding complications.

A P2Y12 inhibitor (i.e., Clopidogrel, Prasugrel, Ticagrelor), must be prescribed post procedure to reduce the risk if thrombosis. Aspirin must be administered concomitantly with the P2Y12 inhibitor, and then continued indefinitely to further reduce the risk of thrombosis.

The optimal duration of DAPT (specifically a P2Y12 platelet inhibitor in addition to aspirin) following DES implantation is unknown, and DES thrombosis may still occur despite continued therapy. It is very important that the patient is compliant with the post-procedural antiplatelet recommendations.

It is strongly advised that the treating physician consider the European Society of Cardiology recommendations for antiplatelet therapy as well as the AHA/ACC/SCAI PCI recommendations for pre- and post-procedure to reduce the risk of thrombosis.

In selected patients at higher risk of bleeding, it may be reasonable to interrupt or discontinue P2Y12 inhibitor therapy after 3 months.

Decisions about duration of DAPT are best made on an individual basis and should integrate clinical judgment, assessment of the benefit/risk ratio, and patient wellbeing.

Premature discontinuation or interruption of prescribed antiplatelet medication could result in a higher risk of stent thrombosis, MI or death.

5.3 USE OF MULTIPLE STENTS

A patient's exposure to drug and polymers is proportional to the number and total length of implanted stents. In the EluNIR™ pivotal clinical trials, total stenting length per subject was limited up to 100 mm.

5.4 USE IN CONJUNCTION WITH OTHER PROCEDURES

The safety and effectiveness of using mechanical atherectomy devices (directional atherectomy catheters, rotational atherectomy catheters) or laser angioplasty catheters in conjunction with EluNIR-PERL™ stent implantation have not been established.

5.5 BRACHYTHERAPY

EluNIR-PERL™ stent safety and effectiveness have not been evaluated in patients with prior target lesion or in-stent restenosis-related brachytherapy.

5.6 USE IN SPECIAL POPULATIONS

5.6.1 PREGNANCY

The EluNIR-PERL™ stent has not been tested in pregnant women or in men intending to father children. Effects on the developing fetus have not been studied. Effective contraception should be initiated before implanting an EluNIR-PERL™ stent and continued for one year after implantation. While there is no contraindication, the risks and reproductive effects are unknown at this time.

5.6.2 LACTATION

A decision should be made whether to discontinue nursing prior to stent implantation, considering the importance of the stent to the mother.

5.6.3 PEDIATRIC USE

The safety and effectiveness of the EluNIR-PERL™ stent in pediatric subjects have not been established. It is not recommended to use the EluNIR-PERL™ stent in a pediatric population.

5.6.4 GERIATRIC USE

The safety and efficacy of EluNIR™ stent in clinical studies has established that the stent can be used in geriatric subject with no upper age limit.

5.7 LESION/VESSEL CHARACTERISTICS

The safety and effectiveness of the EluNIR-PERL™ stent have not been established for subject populations with the following clinical settings:

- Occlusive thrombus and/or a thrombus requiring thrombectomy in a target vessel
- Coronary artery reference vessel diameters <2.25 mm or >4.25 mm
- Unprotected left main lesions ≥30%, or planned unprotected left main intervention
- Ostial LAD or LCX lesions (stenting of any diseased segment within 5 mm of the unprotected left main coronary artery)
- Planned stenting of lesions in more than two (2) major coronary arteries (i.e., two of LAD, LCX, RCA) and their respective branches (the Ramus Intermedius is defined as a branch of the LCX)
- Planned stenting of more than two lesions per vessel (two lesions separated by less than 10 mm and which can be covered by a single stent are considered as one lesion)
- Bifurcation lesions with planned dual stent implantation
- Recent acute myocardial infarction (STEMI)
- Stenting of lesions due to DES restenosis

5.8 DRUG INTERACTIONS

See section 6.2 - Drug Information, Drug Interactions.

Several drugs are known to affect ridaforolimus metabolism, and other drug interactions may also occur. Ridaforolimus is known to be a substrate for both cytochrome P4503A4 (CYP3A4) and P-glycoprotein (PgP). Ridaforolimus absorption and subsequent elimination may be influenced by drugs that affect these pathways. Formal drug interaction studies have not been performed with the EluNIR-PERL™ stent (see section 6.2 - Drug Information, Drug Interactions) because of limited systemic exposure to ridaforolimus eluted from the EluNIR-PERL™ stent. Therefore, due consideration should be given to the potential for both systemic and local drug interactions in the vessel wall, when deciding to place the EluNIR-PERL™ stent in a patient taking a drug with known interaction with ridaforolimus, or when deciding to initiate therapy with such a drug in a patient who has recently received a EluNIR-PERL™ stent.

5.9 MAGNETIC RESONANCE IMAGING (MRI) SAFETY INFORMATION

Non-clinical testing has demonstrated that the EluNIR-PERL™ stent is MR Conditional for single and overlapping stents up to 120 mm. A patient with an implant from this family can be scanned safely in an MR system meeting the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
- Maximum spatial gradient magnetic field of 3,000-gauss/cm (30-T/m) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg in the Normal Operating Mode

Under the scanning conditions defined, the EluNIR-PERL™ stent is expected to produce a maximum temperature rise of 3°C after 15 minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the image artifact caused by the EluNIR-PERL™ stent extends approximately 8 mm from this EluNIR-PERL™ stent when imaged with a gradient echo pulse sequence and a 3-Tesla MR system. The artifact obscures the device lumen..

5.10 STENT HANDLING

- **Each stent is for single use only.** Do not re-sterilize or reuse this device. Note the “Use by” (expiration) date on the product label.
- **The foil pouch is not a sterile barrier.** The inner header bag pouch within the foil pouch is the sterile barrier. **Only the contents of the inner sterilization pouch should be considered sterile. The outside surface of the inner sterilization pouch is NOT sterile.**
- **Do not remove the stent from the delivery system.** Removal may damage the stent and/or lead to stent embolization. These components are intended to perform together as a system.
- The delivery system should not be used in conjunction with other stents.
- Special care must be taken not to handle or disrupt the stent on the balloon, particularly during delivery system removal from packaging, placement over the guide wire and advancement through the rotating hemostatic valve adapter and guiding catheter hub. Manipulation, e.g. rolling the mounted stent with your

fingers, may loosen the stent from the delivery system balloon and cause dislodgment.

- When loading the catheter on the guidewire, provide adequate support to shaft segments.
- **Do not manipulate, touch, or handle the stent with your fingers,** which may cause coating damage, contamination, or stent dislodgement from the delivery balloon.
- Use only the appropriate balloon inflation media (see section 10.4.3 - *Operator's Manual, Delivery System Preparation*). Do not use air or any gaseous medium to inflate the balloon as this may cause uneven expansion and difficulty in stent deployment.

5.11 STENT PLACEMENT

5.11.1 STENT PREPARATION

- **Do not prepare or pre-inflate the delivery system prior to stent deployment other than as directed.** Use the balloon purging technique described in section 10.4.3 - *Operator's Manual, Delivery System Preparation*.
- **While introducing the delivery system into the vessel, do not induce negative pressure on the delivery system.** This may cause dislodgement of the stent from the balloon.
- Use guiding catheters with lumen sizes suitable to accommodate the introduction of the stent delivery system (see section 1 - *Product Description*).

5.11.2 STENT IMPLANTATION

- The decision to pre-dilate the lesion with an appropriate sized balloon should be carefully based on patient and lesion characteristics. The EluNIR™ pivotal clinical trials demonstrated that in a real-world setting, direct stenting with EluNIR-PERL™ stents in single-lesion treated patients who did not have a staged procedure is safe.
- Do not expand the stent if it is not properly positioned in the vessel (see section 5.12 - *Precautions, Stent System Removal*).
- Implanting a stent may lead to a dissection of the vessel distal and/or proximal to the stented portion and may cause acute closure of the vessel, requiring additional intervention (Coronary artery bypass grafting [CABG], further dilatation, placement of additional stents, or other).
- When treating multiple lesions, consider stenting the distal lesion first, followed by stenting of the proximal lesion. Stenting in this order obviates the need to cross the proximal stent when placing the distal stent and reduces the chances of dislodging the proximal stent.
- Additional expansion of a deployed stent may cause a flow limiting dissection. This may be treated by implantation of another stent. When multiple stents are implanted, the ends should overlap slightly.
- Stent placement may compromise side branch patency.
- **Do not exceed Rated Burst Pressure (RBP) as indicated on product label.**
- Balloon pressures should be monitored during inflation. Applying pressures higher than those specified on the product label may result in a ruptured balloon with possible arterial damage and dissection. The stent inner diameter should approximate 1.1 times the reference diameter of the vessel.
- An unexpanded stent may be retracted into the guiding catheter one time only. An unexpanded stent should not be reintroduced into the artery once it has been pulled back into the guiding catheter. Subsequent movement in and out through the distal end of the guiding catheter should not be performed, as the stent may be damaged when retracting the undeployed stent back into the guiding catheter.
- Should **resistance** be felt **at any time** during either lesion access or withdrawal of the pre-stent implantation, the system should be removed per instructions in section 5.12 - *Precautions, Stent System Removal*.
- Stent retrieval methods (e.g., using additional wires, snares, and/or forceps) may result in additional trauma to the coronary vasculature and/or the vascular access site. Complications may include bleeding, hematoma, pseudoaneurysm, or vessel perforation.
- Although the stent delivery system balloon is strong enough to expand the stent without rupture, a circumferential balloon tear distal to the stent and prior to complete stent expansion, could cause the balloon to become tethered to the stent, requiring surgical removal. In case of balloon rupture, it should be withdrawn and, if necessary, a new dilatation catheter exchanged over the guidewire to complete the expansion of the stent.
- Ensure the stented area covers the entire lesion/dissection site and that no gaps exist between stents.
- Do not torque the catheter more than one (1) full turn.
- If reinserting the catheter, flush the guidewire lumen before reinsertion

5.12 STENT SYSTEM REMOVAL

5.12.1 STENT DELIVERY SYSTEM REMOVAL PRIOR TO STENT DEPLOYMENT

If removal of a stent system is required prior to stent deployment, ensure that the guiding catheter is coaxially positioned relative to the stent delivery system, and cautiously withdraw the stent delivery system into the guiding catheter. Should unusual resistance be felt at any time when withdrawing the stent towards the guiding catheter, the stent delivery system and the guiding catheter should be removed as a single unit. This should be done under direct visualization with fluoroscopy.

5.12.2 WITHDRAWAL OF THE STENT DELIVERY SYSTEM FROM THE DEPLOYED STENT

1. Deflate the balloon by pulling negative on the inflation device. Larger and longer balloons will take more time (up to 30 seconds) to deflate than smaller and shorter balloons. Confirm balloon deflation under fluoroscopy and wait 10-15 seconds longer.
2. Position the inflation device to “negative” or “neutral” pressure.
3. Stabilize guiding catheter position just outside coronary ostium and anchor in place. Maintain guide wire placement across stent segment.
4. Gently remove the stent delivery system with slow and steady pressure.
5. Confirm adequate sealing of the hemostatic valve.

If during withdrawal of the delivery system resistance is encountered, use the following steps to improve balloon rewrap:

- Re-inflate the balloon up to nominal pressure.
- Repeat steps 1 through 5 above.

Failure to follow these steps and/or applying excessive force to the delivery system can potentially result in loss of or damage to the stent and/or delivery system components.

If it is necessary to retain guidewire position for subsequent artery/lesion access, leave the guidewire in place and remove all other system components.

Stent retrieval methods (i.e., additional wires, snares, and/or forceps) may result in additional trauma to the coronary vasculature and/or the vascular access site. Complications may include, but are not limited to bleeding, hematoma, pseudoaneurysm, or vessel perforation.

5.13 POST PROCEDURE

- When crossing a newly deployed stent with an intravascular ultrasound (IVUS) catheter, a coronary guidewire, an optical coherence tomography (OCT) catheter, a balloon catheter or delivery system, exercise care to avoid disrupting the stent placement, apposition, and geometry.
- Post-dilatation: All efforts should be made to assure that the stent is not under-dilated. If the deployed stent is not fully apposed to the vessel wall, the stent may be expanded further with a larger diameter balloon that is slightly shorter (about 2 mm) than the stent. The post-dilatation can be done using a low-profile, high pressure, non-compliant balloon catheter. The balloon should not extend outside of the stented region. Do not use the stent delivery balloon for post-dilatation.
- Antiplatelet therapy should be administered post-procedure [see section 5.2 - *Precautions, Antiplatelet Regimen*]. Patients who require early discontinuation of antiplatelet therapy (e.g., secondary to active bleeding) should be monitored carefully for cardiac events. At the discretion of the patient's treating physician, the antiplatelet therapy should be restarted as soon as possible.
- If the patient requires imaging, see section 5.9 - *Magnetic Resonance Imaging (MRI) Safety Information*.

6 DRUG INFORMATION

6.1 MECHANISM OF ACTION

The mechanism by which the EluNIR-PERL™ stent inhibits neointimal growth as seen in pre-clinical and clinical studies has not been established. At the cellular level, ridaforolimus inhibits growth factor-stimulated cell proliferation. At the molecular level, ridaforolimus forms a complex with the cytoplasmic protein FKBP-12 (FK 506 Binding Protein). This complex binds to and interferes with FKBP-12 Rapamycin Associated Protein (FRAP), also known as mammalian target of rapamycin (mTOR), leading to inhibition of cell metabolism, growth, and proliferation by arresting the cell cycle at the late G1 stage.

6.2 DRUG INTERACTIONS

Ridaforolimus is extensively metabolized by the cytochrome P4503A4 (CYP3A4) in the liver, and is a substrate for the counter transporter P-glycoprotein (PgP). Therefore, absorption and subsequent elimination of ridaforolimus may be influenced by drugs that also affect CYP3A4 and PgP pathways. Formal drug interaction studies have not been performed with the EluNIR-PERL™ stent because of limited systemic exposure to ridaforolimus eluted from EluNIR-PERL™. However, consideration should be given to the potential for both systemic and local drug interactions in the vessel wall when deciding to place the EluNIR-PERL™ stent in a subject taking a drug with known interaction with ridaforolimus.

Ridaforolimus, when prescribed as an oral medication, may interact with CYP3A4/PgP inhibitors and/or CYP3A4/PgP inducers. Medications that are strong inhibitors of CYP3A4 or PgP might reduce ridaforolimus metabolism in vivo. Hence, co-administration of strong inhibitors of CYP3A4 or PgP may increase the blood concentrations of ridaforolimus. Medications that are strong inducers of CYP3A4 or PgP might increase ridaforolimus metabolism in vivo resulting in decreased blood concentrations of ridaforolimus.

6.3 MUTAGENESIS, CARCINOGENICITY AND REPRODUCTIVE TOXICITY

The EluNIR-PERL™ stent uses an identical coating (including drug and polymers) as EluNIR™. Therefore, the data from Mutagenesis, Carcinogenicity and Reproductive Toxicity studies of the EluNIR™ stent are considered to be representative of the EluNIR-PERL™ stent.

6.3.1 MUTAGENESIS

The mutagenic potential of the EluNIR™ stent was evaluated in three separate assays, in bacteria and in mammalian cells. The potential to reverse mutations in *S. typhimurium* and *E. coli* strains was assessed in vitro in the standard Ames assay. Mouse lymphoma L5178Y/Tk⁺ cells were also used to evaluate, in vitro, the potential to induce forward mutations at the thymidine kinase (Tk) locus. Lastly, the potential clastogenic effects of the EluNIR™ stent were determined in vivo using the Mouse Micronucleus Test. No evidence of mutagenic potential was observed in these studies.

6.3.2 CARCINOGENICITY

Formal carcinogenicity testing was not conducted on the EluNIR™ Stent. The carcinogenic potential of the EluNIR™ stent is minimal based on the limited period of ridaforolimus release, on the types and quantities of materials present, and on the favorable outcomes and results of the mutagenesis testing.

6.3.3 REPRODUCTIVE TOXICITY

Formal reproductive toxicity testing was not conducted on the EluNIR™ Stent. The reproductive potential of the EluNIR™ stent is minimal based on the limited period of ridaforolimus release, on the types and quantities of materials present.

7 POTENTIAL ADVERSE EVENTS/COMPLICATIONS

Below is a list of the potential adverse effects (e.g., complications) associated with the use of the EluNIR-PERL™ Ridaforolimus Eluting Coronary Stent System.

Adverse events (in alphabetical order) which may be associated with percutaneous coronary and treatment procedures, where coronary stents are used in native coronary arteries include, but are not limited to:

- Access site complications*
- Acute myocardial infarction
- Allergic reaction or hypersensitivity
- to stent components or contrast media
- Aneurysm
- Angina pectoris
- Anxiety
- Bleeding complications which may require transfusions or surgical repair
- Need for CABG-emergent or non-emergent Cardiac arrhythmias
- Cardiac failure
- Cardiac tamponade
- Cardiac shock
- Coronary artery complications**
- Death
- Delayed endothelialization
- Distal emboli
- Endocarditis
- Failure to deliver the stent to the intended site
- Fever or pyrogenic reactions
- Hypertension
- Hypotension
- Infections
- Myocardial ischemia
- Nausea and vomiting
- Palpitations
- Perforation of the heart or great vessels
- Pericardial effusion
- Pulmonary failure
- Renal failure
- Stent compression
- Stent misplacement/migration/embolization
- Stent thrombosis
- Stroke/cerebrovascular accident (CVA) / transient ischemic attack (TIA)
- Vasovagal reaction
- Ventricular fibrillation
- Volume overload

* Includes arteriovenous fistula, hematoma, infection, nerve injury, pain, peripheral ischemia, phlebitis, pseudoaneurysm

** Includes abrupt closure, dissection, embolism, injury, perforation, plaque rupture/shift, restenosis, rupture, spasm, thrombosis, total occlusion

Patient exposure to ridaforolimus is directly related to the total stent length implanted. The adverse events that have been associated with the intravenous injection of ridaforolimus in humans include but are not limited to those listed below. These events are based on experience with the drug in phase I/II studies conducted by Merck Sharp & Dohme Corp. and Ariad Pharmaceuticals Inc. The systemic exposure in representative phase I studies was in concentrations that are 150 times greater than foreseeable with the EluNIR-PERL™ stent.

Potential adverse events (AEs) and adverse drug events (ADEs) for systemic exposure of ridaforolimus:

- Anemia
- Anorexia
- Alopecia
- Aspartate Aminotransferase increased
- Blood Creatine phosphokinase
- Blood Alkaline Phosphatase increased
- Constipation
- Dehydration
- Diarrhea
- Dysgeusia
- Dermatitis acneiform
- Febrile neutropenia
- Fatigue
- Hyperglycemia
- Hypertriglyceridemia
- Hypokalaemia
- Hypercholesterolaemia
- Hypophosphataemia
- Leukopenia
- Mucosal inflammation
- Nausea
- Nail disorder
- Pneumonia
- Pneumonitis
- Pyrexia
- Pruritus
- Paraesthesia
- Renal failure acute
- Rash
- Stomatitis
- Thrombocytopenia
- Vomiting
- Weight decrease

Up to this revision date, no drug related adverse events have been reported. The manufacturer will continue to monitor at post-market stage for any new/unknown risk or adverse events related to drug.

8 INDIVIDUALIZATION OF TREATMENT

The risks and benefits should be considered for each patient before use of the EluNIR-PERL™. Patient selection factors to be assessed should include a judgment regarding risk of long-term antiplatelet therapy. Stenting is generally avoided in those patients at heightened risk of bleeding (e.g., patients with recently active gastritis or peptic ulcer disease) in which anticoagulation therapy would be contraindicated.

Antiplatelet drugs should be used in combination with the EluNIR-PERL™ stent. See also section 5.2 - *Precautions, Antiplatelet Regimen*, section 5.6 - *Precautions, Use in Special Populations*, and section 5.7 - *Precautions, Lesion / Vessel Characteristics*.

Premorbid conditions that increase the risk of poor initial results or the risks of emergency referral for bypass surgery (diabetes mellitus, renal failure, and severe obesity) should be reviewed.

8.1 PATIENT COUNSELING AND PATIENT INFORMATION

Physicians should consider the following in counseling patients about this product:

- Discuss the risks associated with stent placement
- Discuss the risks associated with an ridaforolimus eluting stent
- Discuss the risks of early discontinuation of the antiplatelet therapy
- Discuss the risks of late stent thrombosis with DES use in higher risk patient subgroups
- Discuss the risk/benefit issues for this particular patient
- Discuss alternation to current life style immediately following the procedure and over the long term

8.2 CLINICAL BENEFITS

The primary clinical benefits include the following:

The EluNIR-PERL™ Ridaforolimus Eluting Coronary Stent System maintains vessel patency post stent implantation and reduces stenoses recurrence, resulting in decreased repeat revascularization (a therapeutic effect of the drug ridaforolimus); improving myocardial perfusion, thus alleviating related anginal symptoms. The measurable patient-relevant outcomes quantitative clinical data from Medinol-sponsored EluNIR clinical trials are shown in **Table 2**.

TABLE 2 - CLINICAL BENEFITS

Patient Relevant outcomes	BIONICS HBR Study (N=315 subjects) (N=374 lesions)	BIONICS & NIREUS Pooled Study (N=1159 subjects) (N=1482 lesions)	BIONICS Small Vessels (N=81 subjects) (N=87 lesions)
SAFETY MEASURES			
Clinically-driven target vessel revascularization (TVR) at 12 months ^a	1.9%	4.7%	0.0%
Clinically-driven target lesion revascularization (TLR) at 12 months ^b	1.0%	3.2%	0.0%
EFFICACY MEASURES			
Lesion success ^c	98.7%	99.9%	100%
Device Success ^d	98.7%	98.2%	98.8%

^a Revascularization in the target vessel associated with positive functional ischemia study or ischemic symptoms and an angiographic minimal lumen diameter stenosis $\geq 50\%$ by QCA, or revascularization of a target vessel with diameter stenosis $\geq 70\%$ by QCA without either angina or a positive functional study.

^b Revascularization at the target lesion associated with positive functional ischemia study or ischemic symptoms and an angiographic minimal lumen diameter stenosis $\geq 50\%$ by quantitative coronary angiography (QCA), or revascularization of a target lesion with diameter stenosis $\geq 70\%$ by QCA without either angina or a positive functional study.

^c The attainment of $<30\%$ residual stenosis by QCA (or $<20\%$ by visual assessment) and thrombolysis in myocardial infarction (TIMI) flow 3 after the procedure, using any percutaneous method.

^d The attainment of $<30\%$ residual stenosis by QCA (or $<20\%$ by visual assessment) and TIMI flow 3 after the procedure, using the assigned device only.

9 HOW SUPPLIED

Sterile - This device is sterilized with ethylene oxide gas and is non-pyrogenic. It is intended for single use only. Do not resterilize. Do not use if the package is opened or damaged. Use prior to the expiration date ("Use By" date).

Contents - One (1) EluNIR-PERL™ Ridaforolimus Eluting Coronary Stent System.

Storage - Do not remove from carton until ready for use. Store in a dry place, away from direct sunlight, at up to 25°C (77°F), with transient excursions allowed up to 40°C (104°F).

10 OPERATOR'S MANUAL

10.1 INSPECTION PRIOR TO USE

Carefully inspect the sterile package before opening and check for damage to the sterile barrier (sterilization pouch).

Do not use if the integrity of the sterile package has been compromised.

- Do not use after the "Use by" date.
- Tear open the foil pouch and remove the inner pouch.

NOTE

The outside of the inner sterilization pouch is NOT sterile.

Open the inner sterilization pouch and pass or drop the product into the sterile field using an aseptic technique.

- Prior to using the EluNIR-PERL™, carefully remove the system from the package and inspect for bends, kinks, and other damage. Verify that the stent does not extend beyond the radiopaque balloon markers. Do not use if any defects are noted. However, do not manipulate, touch, or handle the stent with your fingers, which may cause coating damage, contamination, or stent dislodgement from the delivery balloon.

NOTE

At any time during use of the EluNIR-PERL™ RX, if the stainless steel proximal shaft has been bent or kinked, do not continue to use the catheter.

10.2 MATERIALS REQUIRED

The following is a list of the materials required for the procedure:

- Appropriate guiding catheter(s). See Table 1 - *EluNIR-PERL™ Ridaforolimus Eluting Coronary Stent System Product Description*.
- 2-3 syringes (10-20 ml)
- 1,000 u/500 ml sterile heparinized normal saline (HepNS)
- 0.014" (0.36 mm) x 175 cm (minimum length) guidewire
- Rotating hemostatic valve
- Contrast diluted 1:1 with heparinized normal saline
- Inflation device
- Pre-deployment dilatation catheter
- Three-way stopcock
- Torque device
- Guidewire introducer
- Appropriate arterial sheath
- Appropriate anticoagulation and antiplatelet drugs

10.3 INCIDENTS REPORT

In case of a serious incident related to the EluNIR-PERL™, please contact Medinol at complainthandling@medinol.com and report to the local health authority

10.4 PREPARATION

10.4.1 PACKAGING REMOVAL

NOTE

The foil pouch is not a sterile barrier. The inner header bag (sterilization pouch) within the foil pouch is the sterile barrier. Only the contents of the inner pouch should be considered sterile. The outside surface of the inner sterilization pouch is NOT sterile.

1. Carefully remove the delivery system from its protective tubing for preparation of the delivery system. When using a Rapid Exchange (RX) system, do not bend or kink the hypotube during removal.
2. Remove the product mandrel and protective stent sheath by grasping the catheter just proximal to the stent (at the proximal balloon bond site), and with the other hand, grasp the stent protector and gently remove distally. If unusual resistance is felt during product mandrel and stent sheath removal, do not use this product and replace with another. Follow product returns procedure for the unused device.

10.4.2 GUIDEWIRE LUMEN FLUSH

Flush the guidewire lumen with HepNS until fluid exits the distal tip.

NOTE

Avoid manipulation of the stent while flushing the guidewire lumen, as this may disrupt the placement of the stent on the balloon.

10.4.3 DELIVERY SYSTEM PREPARATION

1. Prepare an inflation device/syringe with diluted contrast medium.
2. Attach an inflation device/syringe to the stopcock; attach it to the inflation port of the product. Do not bend the product hypotube when connecting to the inflation device/syringe.
3. With the tip down, orient the delivery system vertically.
4. Open the stopcock to the delivery system; pull negative for 30 seconds and then release to neutral for contrast fill.
5. Close the stopcock to the delivery system; purge the inflation device/syringe of all air.
6. Repeat steps 3 through 5 until all air is expelled. If bubbles persist, do not use the product.
7. If a syringe was used, attach a prepared inflation device to the stopcock.
8. Open the stopcock to the delivery system.
9. Leave on neutral.

NOTE

While introducing the delivery system into the vessel, do not induce negative pressure on the delivery system. This may cause dislodgement of the stent from the balloon.

NOTE

If air is seen in the shaft, repeat section 10.4.3 - Operator's Manual, Delivery System Preparation, steps 3 through 5, to prevent uneven stent expansion.

10.5 DELIVERY PROCEDURE

1. Prepare the vascular access site according to standard practice.
2. The decision to pre-dilate the lesion with an appropriate sized balloon should be based on patient and lesion characteristics. **If pre-dilatation is performed**, limit the longitudinal length of pre-dilatation by the PTCA balloon to avoid creating a region of vessel injury that is outside the boundaries of the EluNIR-PERL™ stent.
3. For long lesions, size the stent to the diameter of the most distal portion of the vessel.

NOTE

If choosing between two stent diameters for tight lesions, choose the smaller diameter stent and inflate.

4. Maintain neutral pressure on the inflation device attached to the delivery system. Open the rotating hemostatic valve as wide as possible.
5. Backload the delivery system onto the proximal portion of the guidewire while maintaining guidewire position across the target lesion.
6. Carefully advance the delivery system into the guiding catheter and over the guidewire to the target lesion. When using a Rapid Exchange (RX) system, be sure to keep the hypotube straight. Ensure guiding catheter stability before advancing the stent system into the coronary artery.

NOTE

If unusual resistance is felt before the stent exits the guiding catheter, do not force passage. Resistance may indicate a problem and the use of excessive force may result in stent damage or dislodgement. Maintain guidewire placement across the lesion and remove the delivery system and guiding catheter as a single unit.

7. Advance the delivery system over the guidewire to the target lesion under direct fluoroscopic visualization. Utilize the radiopaque balloon markers to position the stent across the lesion. Perform angiography to confirm stent position. If the position of the stent is not optimal, the stent should be carefully repositioned or removed [see section 5.12 - *Precautions, Stent System Removal*]. The balloon markers indicate both the stent edges and the balloon shoulders. Expansion of the stent should not be undertaken if the stent is not properly positioned in the target lesion.

NOTE

If removal of a stent system is required prior to deployment, ensure that the guiding catheter is coaxially positioned relative to the stent delivery system, and cautiously withdraw the stent delivery system into the guiding catheter. Should unusual resistance be felt at any time when withdrawing the stent towards the guiding catheter, the stent delivery system and the guiding catheter should be removed as a single unit. This should be done under direct visualization with fluoroscopy.

8. Tighten the rotating hemostatic valve. The stent is now ready to be deployed.

10.6 DEPLOYMENT PROCEDURE

1. Prior to deployment, reconfirm the correct position of the stent relative to the target lesion using the radiopaque balloon markers.
2. Deploy the stent slowly by pressurizing the delivery system in 2 atm increments, every 5 seconds, until stent is completely expanded. Fully expand the stent by inflating to nominal pressure at a minimum. Accepted practice generally targets an initial deployment pressure that would achieve a stent inner diameter ratio of about 1.1 times the reference vessel diameter.
3. For long lesions, size the stent to the diameter of the most distal portion of the vessel and expand the stent to nominal pressure at minimum. Maintain pressure for 30 seconds. If necessary, the delivery system can be re-pressurized or further pressurized to assure complete apposition of the stent to the artery wall.
4. Maintain pressure for 30 seconds for full expansion of the stent. Fluoroscopic visualization during stent expansion should be used in order to properly judge the optimum stent diameter as compared to the proximal and distal native coronary artery diameters [reference vessel diameters]. Optimal stent expansion and proper apposition requires that the stent be in full contact with the arterial wall.

NOTE

See section 10.7 - *Removal Procedure* for instructions on withdrawal of the stent delivery system.

5. If necessary, the delivery system can be re-pressurized or further pressurized to assure complete apposition of the stent to the artery wall.

NOTE

Do not exceed the labeled rated burst pressure (RBP) of 18atm (1824kPa).

6. Fully cover the entire lesion and balloon treated area (including dissections) with the EluNIR-PERL™ stent, allowing for adequate stent coverage into healthy tissue proximal and distal to the lesion.
7. Deflate the balloon by pulling negative on the inflation device for 30 seconds. Confirm complete balloon deflation before attempting to move the delivery system. If unusual resistance is felt during stent delivery system withdrawal, pay particular attention to guiding catheter position.

NOTE

See section 10.7 - *Removal Procedure* for instructions on withdrawal of the stent delivery system.

8. Confirm stent position and deployment using standard angiographic techniques. For optimal results, the entire stenosed arterial segment should be covered by the stent. Fluoroscopic visualization during stent expansion should be used in order to properly judge the optimum expanded stent diameter as compared to the proximal and distal coronary artery diameter(s). Optimal expansion requires that the stent be in full contact with the artery wall. Stent wall contact should be verified through routine angiography or intravascular ultrasound (IVUS).
9. If the deployed stent size is still inadequate with respect to reference vessel diameter, a larger balloon may be used to further expand the stent. If the initial angiographic appearance is suboptimal, the stent may be further expanded using a low profile, high pressure, non-compliant balloon dilatation catheter. If this is required, the stented segment should be carefully re-crossed with a prolapsed guidewire to avoid disrupting the stent geometry. **Deployed stents should not be left under-dilated.**

CAUTION

Do not dilate the stent beyond the following limits:

NOMINAL STENT DIAMETER	DILATATION LIMITS
2.25 mm	3.00 mm
2.5 mm	3.00 mm
2.75 mm	3.75 mm
3.5 mm	3.75 mm
3.0 mm	4.75 mm
4.0 mm	4.75 mm

10.7 REMOVAL PROCEDURE

10.7.1 WITHDRAWAL OF THE STENT DELIVERY CATHETER FROM THE DEPLOYED STENT:

1. Deflate the balloon by pulling negative on the inflation device. Larger and longer balloons will take more time (up to 30 seconds) to deflate than smaller and shorter balloons. Confirm balloon deflation under fluoroscopy and wait 10-15 seconds longer.
2. Position inflation device on "negative" or "neutral" pressure.
3. Stabilize the guiding catheter position just outside coronary ostium and anchor in place. Maintain guidewire placement across the stent segment.
4. Gently remove the stent delivery system with slow and steady pressure.
5. Confirm adequate sealing of the hemostatic valve..

If resistance is encountered during withdrawal of the stent delivery catheter, use the following steps to improve balloon rewrap:

6. Re-inflate the balloon up to nominal pressure.
7. Repeat steps 1 through 5 above.

10.7.2 POST STENT DELIVERY SYSTEM WITHDRAWAL - STENT DEPLOYMENT CONFIRMATION:

1. Confirm stent position and deployment using standard angiographic techniques. For optimal results, the entire stenosed arterial segment should be covered by the stent. Fluoroscopic visualization during stent expansion should be used in order to properly judge the optimum expanded stent diameter as compared to the proximal and distal coronary artery diameter(s). Optimal expansion requires that the stent be in full contact with the artery wall. Stent wall contact should be verified through routine angiography or intravascular ultrasound (IVUS).
2. If more than one EluNIR-PERL™ stent is needed to cover the lesion and balloon treated area, it is suggested that, to avoid the potential for gap restenosis, the stents be adequately overlapped.
3. To ensure that there are no gaps between stents, the balloon marker bands of the second EluNIR-PERL™ stent should be positioned inside the deployed stent prior to expansion.
4. Reconfirm stent position and angiographic results to assess stented area. Repeat inflations until optimal stent deployment is achieved. If post-dilatation is necessary, ensure that the final stent diameter matches the reference vessel diameter. **Assure that the stent wall is in contact with the artery wall.**

10.8 POST-DEPLOYMENT DILATATION OF STENT SEGMENTS

1. All efforts should be made to assure that the stent is not under-dilated.
2. If the deployed stent size is still inadequate with respect to the vessel diameter, or if full contact with the vessel wall is not achieved, a larger balloon may be used to expand the stent further. This is done by using a low profile, high pressure, and noncompliant balloon catheter. If this is required, the stented segment should be re-crossed carefully with a prolapsed guidewire to avoid dislodging the stent. The balloon should be centered within the stent and should not extend outside of the stented region).

CAUTION

Do not dilate the stent beyond the dilatation limits described in section 10.6-9.

10.9 PATIENT IMPLANT CARD & PATIENT INFORMATION GUIDE

The physician or hospital shall complete the patient, hospital, physician and implant information on the patient implant card. The patient shall be provided with the entire patient implant card and patient information guide before being discharged.

11 REUSE PRECAUTION STATEMENT

Do not use if sterile barrier is damaged.

If damage is found, call your Medinol representative.

For single patient use only.

Do not reuse, reprocess, or resterilize since this may compromise device performance, cause contamination and/or result in severe patient injury.

Safe Disposal.

After use, dispose of the product and packaging in accordance with hospital, administrative and/or local government policy.

12 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

A summary of the safety and clinical performance (SSCP) for this device is available at <https://ec.europa.eu/tools/eudamed>. This is the SSCP location after the launch of the European Database on Medical Devices (EUDAMED).

13 DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Medinol product(s) described in this publication. Under no circumstances shall Medinol be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind Medinol to any representation or warranty except as specifically set forth herein.

Descriptions or specifications in Medinol printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Medinol will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

EluNIR-PERL™ is a trademark of Medinol Ltd.

Medinol's products are covered by one or more pending and issued European and United States patents.

EluNIR-PERL™

Ridaforolimus Eluting Coronary Stent System

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EluNIR-PERL™ je ochranná známka společnosti Medinol Ltd.